

TURIOS

Our most compact






*Scan this QR code
and discover Turios
in the video.*



■ The most compact

Turios is our most compact powered wheelchair. Its small size ensures easy manoeuvrability, making it highly suitable for indoor use. Taking public transport or using a lift is no problem at all with this powered wheelchair.

When used outdoors, the performance of Turios is excellent. With its large wheels, Turios is well suited to overcome small obstacles while keeping optimal stability at all times. Turios, your best trip companion, for indoor and outdoor use!



■ Customize your Turios

Configure your Turios according to your lifestyle and taste! Whether you want to stay most of your time indoors or venture outdoors, your Turios adapts to any situation. The accentuated wheel rim will always follow the color of the Turios frame that you selected during your order.

Front wheels



8" (200 mm)



9" (230 mm)



10" (260 mm)

Rear wheels



■ Configure your Turios

Comfort is very important in a wheelchair which is why you will have several choices of seats and backrests when ordering your Turios. For better positioning and a more stable posture, we can recommend a preshaped backrest or cushion either with PU Soft Touch finishing or stylish artificial leather.



*Turios C5 black - Headrest L55
LED lights SE19 - Comfort leg support BZ8*

Preshaped backrest and seat L34 L35 PU Soft Touch



*Preshaped backrest
L34 3D Net*



*Preshaped backrest
L34 Artificial leather*



*Preshaped backrest
L34 Nylon Plus*



*Preshaped seat
L35 3D Net Nylon Plus*



*Preshaped seat
L35 PU Soft Touch*




*Preshaped seat
L35 Soft artificial leather*

For comfortable seating, Turios is equipped with a manually adjustable seat and backrest or with an optional electrical adjustable backrest and seat providing both comfort and ease of use.



Tension-adjustable backrest - L44



■ Travel plans? Bring your Turios along!



The Turios' batteries are allowed to be transported by air making it easy to bring your Turios with you to your favorite destination.



Turios has been crash tested according to **ISO 7176-19** and can be used as a forward-facing seat in motorized vehicles by children and adults with a body mass higher or equal to 22 kg.

■ Colours


Turios is available in more than 7 different colours. Your chosen colour will be applied on the lower frame and also on the rims of the rear wheels.




C4
red



C5
black



C76
metallic white



C77
metallic black




C79
metallic bright blue



C85
green



C87
orange



Other colors
on request



SE09 - Attendant controller



SE06 - USB connection



B31 - Crutch holder



BZE - Electrical articulated legrests



SE19 - Lights



SE27 - Top Case

LINX™


LINX, The next generation wheelchair control system

- A seamless and expandable system
- Modular components – mix and match modules
- State-of-the-art technology
- Revolutionary approach to programming



A complete list of all possible options is available on the website www.vermeiren.com

Technical data

■ Total width	580 - 640 mm	■ Batteries	2x12V/38Ah/150W or 2x12V/50Ah/450W
■ Total length	1060 - 1070 mm	■ Speed	6 km/h or 10 km/h
■ Seat height	450 - 550 mm	■ Driving range	>25 km
■ Seat width	390 - 450 - 510 mm	■ Weight	105,8 kg
■ Seat depth	380 - 510 mm	■ Maximum user weight	130 kg
■ Backrest height	550 or 585 mm	■ Additional info	

Measurements in mm, weight in kg - Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

TURIOS

Our most compact

- compact design
- indoor and outdoor use
- powerful rear-wheel drive
- outstanding driving comfort

distributed by:

VERMEIREN
we care for you

