


Material				Material			
No.	Plan	Article	Quantity	No.	Plan	Article	Quantity
1	YP09307DP	26435...	1	34	XP00071	2990040	2
2	YP00372	18117...	1	35	Screw D4x12	1990816	2
3	Y004067	15072...	1	36	X004865	1003901	4
4	Y004072	15073...	2	37	---	1994587	2
5	YP09302	2643404	2	38	X002335	1993411	1
6	---	1611704	1	39	M5x10 DIN 916/45H	1001485	1
7	Front Wheel	1916560	2	40	L=980mm	1003279	2
8	M12x60	1993419	2	41	D6,4 DIN125A	1001564	2
9	Forks	1991447	2	42	XP00028	2990012	2
10	M10x80 DIN931	1001457	2	43*	---	PP00381	2
11	M10 DIN 985	1001683	2	44*	LOGO VERMEIREN	T001002	1
12	M6x75 DIN 912/8.8	1993413	2	45	M6x35 DIN931	2995333	2
13	X002306	1502304	2	46*	---	1003558	1
14	XP01452	2706000	2	47*	---	W620150	2
15	X002448	1517200	2	48*	---	W620151	2
16	D6,4 DIN 125 A	1001564	4	49	XP00108	GP00604	2
17	---	1003558	2	50	M8x76 DIN912/8.8	1998993	2
18	M6 DIN985	1001680	10	51	M8	1998994	2
19	M6x40	1006111	2	52	Rear Wheel	1914606	2
20	M6x60 DIN 912/8.8	1005957	4				
21	M6 DIN125A	1000973	4				
22	XP00099	2990067	4				
23	D10,5 DIN 125 A	1001558	6				
24	M6x16	1008459	2				
25	X002305	1001545	4				
26	X002305	1993120	2				
27	D5x25	H016007	2				
28	X005513	1992600	2				
29**	YP05288	2993604	2				
30	---	1993417	1				
31	Sterknop M6	2995334	2				
32	M12 DIN 127 A	1001572	2				
33	Handle and Line L=960mm	1909113	1				

* - Not depicted

** - For Switz only

*** - Option FX + carton U000254/

EXPLODED VIEW 275				Plan No.
2/2				E000042/O
Scale:	Drafted	Verified	Confirmed	
Date	15/04/2003	06/12/2022	06/12/2022	
Init.	BS	PM	AZ	